

Disposable Cup Electrodes

The disposable cup electrodes are one of the latest additions to our product programme.

The diameter of the cups is 10 mm and they are made of silver-chlorided plastic. The PVC wires connects into a safety DIN 42 802 connector, are multi colored (5) and available in the following lengths: 1, 1.5 or 2 meters.

For single use only.
Packing: 25 pieces per package.

Wire Length	Part Number
1 m	C52-438
1.5 m	C52-638
2 m	C52-838



Reusable Cup Electrodes

Our reusable cup electrodes have been designed so they can be used repeatedly by medical technicians.

The precision stamped cups are available in two diameters: the pediatric, 6 mm or the adult, 10 mm size.

They are made out of 99.9% pure silver. The silver cups are silver-chlorided, covered with a heavy gold-plate, or kept silver.

The flexible autoclavable wires connects into a safety DIN 42 802 connector, multi color (6) and are available in the following lengths: 1, 1.5 or 2 meters.

Packing: 12 pieces per package.

Material	Size	Wire Length	Part Number
Silver (Ag)	Adult	1 m	C12-434
Silver (Ag)	Adult	1.5 m	C12-634
Silver (Ag)	Adult	2 m	C12-834
Silver (Ag)	Pediatric	1 m	C14-434
Silver (Ag)	Pediatric	1.5 m	C14-634
Silver (Ag)	Pediatric	2 m	C14-834
Silver-Silver Chloride (Ag-AgCl)	Adult	1 m	C22-434
Silver-Silver Chloride (Ag-AgCl)	Adult	1.5 m	C22-634
Silver-Silver Chloride (Ag-AgCl)	Adult	2 m	C22-834
Silver-Silver Chloride (Ag-AgCl)	Pediatric	1 m	C24-434
Silver-Silver Chloride (Ag-AgCl)	Pediatric	1.5 m	C24-634
Silver-Silver Chloride (Ag-AgCl)	Pediatric	2 m	C24-834
Gold-plated (Au)	Adult	1 m	C32-434
Gold-plated (Au)	Adult	1.5 m	C32-634
Gold-plated (Au)	Adult	2 m	C32-834
Gold-plated (Au)	Pediatric	1 m	C34-434
Gold-plated (Au)	Pediatric	1.5 m	C34-634
Gold-plated (Au)	Pediatric	2 m	C34-834



Moulded safety DIN 42 802 connector and cup ends

Silver (Ag)	Adult	1.5 m	C12-634/01
Gold-plated (Au)	Adult	1.5 m	C32-634/01

