

Disposable Subdermal Corkscrew Needle Electrodes

Stainless Steel corkscrew needle 0.6 mm (23G) in diameter.

The multi color (6) PVC wires connect into a safety DIN 42 802 connector and are available in a 1.2 meter length.

4 x 6 colors (Red, Green, Black, White, Yellow, Blue), 24 pieces per box.
Sterilization: EtO. For single use only.

Wire Length	Part Number
1.2 m	S50715-001



Disposable Subdermal Single Needle Electrodes

Standard single Stainless steel needle, 13 mm in length and 0.4 mm (27G) in diameter.

The multi color (5) PVC wires connect into a safety DIN 42 802 connector and are available in a 1, 1.5, 2 and 2.5 meter length.

4 x 6 colors (Red, Green, Black, White, Yellow, Blue), 24 pieces per box.
Sterilization: EtO. For single use only.

Wire Length	Part Number
1 m	S50716-001
1.5 m	S50716-002
2 m	S50716-003
2.5 m	S50716-004



5 x 5 colors (Red, Green, White, Yellow, Blue), 25 pieces per box.
Sterilization: EtO. For single use only.

Wire Length	Part Number
1 m	S43-438
1.5 m	S43-638
2 m	S43-838
2.5 m	S43-938

Disposable Subdermal Curved Needle Electrodes

The hooked needle is especially designed to better stay in place during usage. The item consists of a stainless steel bent needle with a 1.5 or 2 meter PVC wire, connects into a safety DIN 42 802 connector.

5 x 5 colors (Red, Green, White, Yellow, Blue), 25 pieces per box.
Sterilization: EtO. For single use only.

Wire Length	Part Number
1.5 m	S45-638
2 m	S45-838

