HugeMed Resuable Video Laryngoscope





The image quality

- 4inch full-angle LCD high resolution display, Resolution ratio: 800*480RGB
- Rotation: up/down:140° left/right:180°



Functional

- one-click button for photography and video recording
- built-in 8G storage card



Model is complete

 8 models of blades are disigned for various population

HugeMed Disposable Video Laryngoscope







Convenient for teaching

- Through the USB data interface, charging, video copy the output
- Extensible external HD interface display



High resolution wide-angle

- 2 megapixels camera
- Field angle ≧60 degrees, and the lens angle is between 50°-60°



Model is complete

 5 models of blades are disigned for various population

Safety

Anti-Fog Lens
Anti-Microbial Handle

Disposable blades

Sterile package
Single use design avoid Cross Infection

Reusable blades

Durable 316 Medical stainless steel IPX8 waterproof grade design Support immersion disinfection

Faster

Lightweight, Portable

Faster Installation and Intubation

Faster Removal of blades or Drop blades

High-Definition

The 3.5 inch anti-reflective display with full view has the high performance visualization capabilities for every intubation

Better View: 2.0 Megapixels Cameras

Easier

Rechargeable li-ion battery, One Shot Button for Recording

Cost Effective

Low cost disposable blades 316 medical stainless reusable blades can be used 1000 times

Technical Specifications

Display	Size: 4inch Full view
	Monitor: LCD
	Resolution Ratio: 800*480RGB
	Aspect Ratio: 4:3
	Video refresh Rata: 30FPS
	Illumination: LED
Data output	Easy to establish and store files
Image output	USB interface
Laryngoscope blade Camera Include	Resolution Ratio 2.0M Pixels
	Field Angle: 66°
	Illuminance ≥800LUX
Handle	Confortable ergonomic short handle design
Battery	Type Rechargeable Lithium Battery
	Long working time ≥120-150 Minutes
	Voltage 3.7V
	Capacity 3200mAh
	Battery life Cycle: More Than 300 times
	Charging times: less than 4–6hrs